Blowing, Filling and Sealing (BFS) System

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Equip global pharmaceutical manufacturing industry with intelligent wings.


无菌塑料包装吹灌封 (BFS) 系统 - 01 -
Blowing, Filling and Sealing (BFS) System

间接式吹灌封 (BFS) 生产线 - 04 -
Batch-type blowing filling and sealing (BFS) production line

连续式吹灌封 (BFS) 生产线 - 07 -
Continuous blowing filling and sealing (BFS) production line

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Founded in 2000, Truking Technology Limited has become one of the major pharmaceutical equipment suppliers in the world. The company specializes in pharmaceutical equipment and the corresponding integrated solutions, and takes the initiative to push forward the R&D of smart pharmaceutical factory. It is a listed A-share company with multiple wholly-owned or controlling subsidiaries, namely Germany based Romaco Group, Truking Watertown, Sichuan Pharmaceutical Design Institute, Truking Feiyun, Truking Ingenuity Biotechnology, Truking Wachine, Truking Intelligent Robot, etc. The company has over 6500 employees around the global, a total asset of RMB 6.9 billion and annual output value of RMB 8 billion.


TRUKING has two operation headquarters respectively in Changsha, China and Germany; 3 big research institutes, namely Changsha Central Technology Research Institute, Europe Technology Research Institute, and Sichuan Pharmaceutical Design Institute. Truking has established National Enterprise Technology Center, National Innovation Center, Post-Doctoral Research Center, Academician Expert Workstation and many more technology and innovation platforms. To better respond to the needs of global customers, Truking has been setting up sales and service centers in many countries including India, Italy, Greece, Russia, Uzbekistan, Turkey, Egypt, Bangladesh, Korea, Vietnam, Thailand, etc. By Jan. 31, 2022, Truking has submitted 4005 patent applications, and been granted with 2427 patents. Beyond China, Truking has applied for 40 PCT international patents, and been granted with 20 patents in USA, Russia, India, Korea, Germany, Indonesia, and many more European countries. Truking also took the lead in formulating 17 national industrial product technical standards. Our products have been exported to more than 40 countries and regions to Asia, Europe, and South America. The international market share is increasing rapidly each year.

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The overall solution BFS TRUKING Technology consists of six core functional modules: BFS equipment, liquid distribution system, sterilization and logistics system, lamp inspection and leak detection system, post packaging system and centralized control and scheduling system. Based on high efficiency and energy saving process equipment system, information system as the core, and automatic logistics as the link to achieve high efficiency, high quality BFS products, safe production.

工艺流程 / BFS Technical Process

挤塑 Extruding	<p>螺杆将粒子塑化连续挤出通过模头形成管胚。管胚内有无菌空气支撑。塑料热融和挤出的温度为 170~220 度，是经过实证的无菌挤出和去热源工艺。</p> <p>The screw extrudes the plasticized particles continuously through the mold end to form the tube embryo. The tube embryo is supported by sterile air. The temperature of plastic melted and extruded is 170~220 degrees</p>
成型 Blowing	<p>模具在挤出工位合模后快速移动到吹瓶 / 灌装工位，A 级风淋装置中对吹瓶、灌装和封口工位进行保护。通过真空和 / 或无菌压缩空气按照模具的形状，把塑料原料成型容器。</p> <p>The mould is quickly moved to the bottle blowing/filling station after the extruding station. The bottle blowing, filling and sealing stations are protected in the A-class air shower device. The plastic material is formed into a container through vacuum and/or sterile compressed air in the shape of mold</p>
灌装 Filling	<p>通过时间压力阀将药液灌装至容器里，并排出无菌空气。药液经过的管路在灌装之前进行在线 CIP 和在线灭菌 SIP，确保药液安全性。</p> <p>Through the time pressure valve, the liquid is filled into the container and sterile air is discharged. The pipeline can be CIP and SIP shall before filling, ensures safety of medicine</p>
封口 Sealing	<p>灌装结束后，模具顶部和开口保持架之间的塑料处于半融化状态，将封口模具合闭，形成容器的顶部，并使瓶子密封。</p> <p>After filling, the plastic between the top of the mold and the open cage is in a semi-melted state. Close the sealing mold to form the top of the container and seal the bottle.</p>
出瓶 Outfeeding	<p>瓶子封口结束后模具打开，已经完成成型、灌装并密封的瓶子被送出机器，进行切边。</p> <p>After sealing, the mould is opened, bottles have been formed, filled and sealed are sent out of the machine for trimming</p>


BFS 产品应用领域 / Application Fields

注射剂 & 口服液
Injection & Oral Liquid

滴眼剂
Eye Drops


吸入制剂
Inhalation

次抛化妆品
Daily use cosmetics

间歇式吹灌封 (BFS) 生产线 Batch-type blowing filling and sealing (BFS) production line

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Blowing, Filling and sealing technology is Blow - Fill - seal , hereinafter referred to as BFS. In the process of continuous operation of one set of equipment, the mold performs reciprocating motion between the extrusion site and the filling site to complete the molding, filling and sealing of the product. All these processes are completed in a single time under aseptic conditions, and the whole cycle is completed within 12-14 seconds.

Blow-fill-seal technology is Blow-fill-seal, referred to as BFS. The technology is to complete the forming, filling and sealing of the product by reciprocating the die in the extrusion station and the filling station in a continuous operation process of one piece of equipment, all of which are completed under aseptic conditions , the whole cycle is completed within 12~14 seconds.

BFS technology has been used for about 40 years in the world, and it has been proved that the contamination rate can reach less than 0.1%, Based on the summary and analysis of the media filling data, the contamination rate of BFS system is 0.001%.

性能优势 / Performance Advantage

- 设备产能：
间歇式 BFS 设备稳定产能 24000pcs/h (2mL), 为国内 BFS 设备最大稳定产能，技术成熟，稳定可靠。
 - 挤出系统：
欧洲研究院设计，精密加工模头，伺服驱动，管胚动态排气补气，保证挤出系统的稳定性。容器壁厚均一性 $\pm 0.05\text{mm}$ 。
 - 成型系统：
模具开合模和封口模具开合模采用伺服液压系统进行 PID 闭环曲线控制，减少设备的冲击和震动，避免模具磨损。伺服摆动油缸驱动高速移模，降低容器敞口暴露的风险。
 - 灌装系统：
时间压力法灌装，采用高频响电磁阀岛控制装量。配置自动废瓶识别技术，混悬液产品均一性偏差 $\pm 0.5\%$ ，残氧量 $< 1\%$
 - 成型模具
模具制造中心精密制造，材料选用进口无毒高导电铜合金加工，保证模具使用寿命，提高产品外观质量。
 - 液压系统：
伺服液压泵站，出口压力 PID 调节，油缸动作闭环曲线控制，手动状态下，液压泵站低速运行，便于设备调试维护以及提高人员的安全保护。
 - 控制系统
采用 SIEMENS 运动控制器和工业电脑控制，实现电子签名和电子批记录，满足数据完整性的验证要求。
- + Equipment capacity
Stable production capacity of 24000 pcs/h (2 mL), is the largest stable production capacity for domestic BFS equipment, mature technology, stable and reliable.
 - + Extruding system:
Designed from Truking (European) R&D Center, precision machining die head, servo drive, tube embryo dynamic exhaust and fill gas to ensure the stability of the extruding system. Container wall thickness uniformity $\pm 0.05\text{mm}$.
 - + Molding system:
A servo hydraulic system is used to control the PID closed-loop curve to reduce the assault and vibration of the equipment and avoid the wear of the mould. Servo swing cylinder drives high-speed mold shifting to reduce the risk of container exposure.
 - + Filling system:
Time pressure method filling, using high frequency solenoid valve control loading. Equipped with automatic identification technology of waste cylinders, the suspension uniformity deviation $\pm 0.5\%$, and the remaining oxygen $< 1\%$
 - + Mold:
Mold manufacturing center adopts the precision manufacturing technology, and selects the imported non-toxic high-conductivity copper alloy for materials processing, to guarantee the service life of the mold and improve the product appearance quality.
 - + Hydraulic systems:
Servo hydraulic pump station, PID adjustment of outlet pressure, closed-loop curve control of cylinder action. In the manual state, the hydraulic pump station runs at low speed, which is convenient for equipment debugging and maintenance and improves the safety protection of personnel.
 - + Control system
Adopt SIEMENS motion controller and industrial computer to control, the action execution is more accurate implementation of electronic signature and electronic batch records to meet the data integrity verification requirements

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Continuous blowing filling and sealing (BFS) production line

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+ Equipment Capacity

The maximum stable production capacity of continuous BFS equipment designed by Truking (European) R&D Center is 33000pcs/h (2ml), and the utilization rate of product materials is up to 83%.

+ Extruding system

For the closed single extrusion parison, the die head is made of Europe-imported alloy materials that are processed by five-axis high-speed machine tool, surface adopts special coating, screw rod adopts servo drive, pipe parison adopts dynamic exhaust and air make-up, so as to maintain the stability of extrusion system and uniformity of container wall thickness.

+ Molding System:

Multi-mold continuous equipment uses multiple sets of molds installed on a continuous rotating chain, and the chain is made of alloy materials and thermally treated to guarantee its service life.

Up-and-down reciprocating molding system (single-mold continuous type equipment) adopts the full servo drive to open, close, lift and lowers the mold, thus ensuring motion control accuracy.

+ Filling System:

The filling pin is installed inside the extrusion die head, protected by sterile air, and used to complete the filling of product in a closed parison. Time pressure method is used for filling, and high-frequency response electromagnetic valve island is adopted to control the filling amount, and the maximum filling accuracy can be up to $\pm 0.01\text{mL}$.

+ Forming Mold

Mold manufacturing center adopts the precision manufacturing technology, and selects the imported non-toxic high-conductivity copper alloy for materials processing, to guarantee the service life of the mold and improve the product appearance quality.

+ Control System

SIEMENS motion controller and industrial computer are used as control unit to realize electronic signature and electronic batch recording, so as to meet the requirements for data integrity verification.

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
Service Concept

400-9988-900

Adhering to concept of “all for customers, for all customers” , and with the mission of “creating value for customers” , TRUKING focuses on building first-class R&D, design, sales, manufacturing, service teams and service management system to offer high-efficient and high-quality products and services for the customers.

◎ Project Planning Support

TRUKING has a series of professional medical and pharmaceutical design institutes, product research institutes, a group of senior scientific research and technical personnel, national-level aseptic pharmaceutical process verification experts, professors, etc., to provide scientific and effective technical planning support for customers' projects.

◎ Service In The Park

As an equipment supplier, TRUKING not only integrates industrial regulations and standards into product design and manufacturing, but also introduces the industrial regulations and standards into Truking industrial park in the form of lectures and exchanges to interact and integrate with the technology and production of pharmaceutical enterprises.

◎ Process Assurance Service

TRUKING launched the overall solution for medical & pharmaceutical equipment realizing the full automation and aseptic production process, intelligent detection and sorting of the drug, as well as computer verification assurance so as to produce assured products for market.

◎ Equipment Management Consultation Service

TRUKING aims to satisfy the maximum benefited demands of the customer and to provide a package of equipment management services to realize minimization investment and maximization production output together with pharmaceutical factories.

◎ Adhere to the philosophy of "customer-centric", TRUKING creates the extraordinary service experience with the spirit of the craftsman, to create more value for customers.

◎ We are always ready to escort the normal operation of the customer's equipment.

- + Onsite equipment installation, commissioning, acceptance
- + Training for equipment process, operation and maintenance
- + Periodic door-to-door maintenance and return visit
- + Performance test and maintenance for machinery parts
- + Timely supply of spare parts
- + Product technology upgrading and secondary development
- + Equipment overhaul

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